



Introducing  
**V30**  
By Viora

In today's competitive business environment, you need to have the most advanced technology to best meet your customers' demands. However the purchase of numerous, different devices to meet these demands can be at great expense, and the technology can become obsolete while you're still paying off your investment.

Viora's V30 Platform is designed to answer those challenges and more, by providing you with the most advanced technologies within one platform, allowing you to expand your practice as you need it and when you're ready.



PLATFORMS



HANDPIECES



TREATMENTS



COMBINATION



SERVICE



**About Viora**

Viora, a leading provider of medical aesthetic solutions, is passionate about helping customers deliver the most effective results to their patients. Viora improves upon treatment experiences by enhancing existing technology to give customers the best non-invasive solutions available. With a presence in over 65 countries, Viora's non-invasive products address a wide range of symptoms, treating today's most commonly requested applications for anti-aging, body & facial contouring, cellulite reduction, acne clearance, stretch mark reduction, hair removal, skin lightening and more.

[www.vioramed.com](http://www.vioramed.com) | [info@vioramed.com](mailto:info@vioramed.com)

**Manufacturer's Address  
Viora Ltd.**

3 Maskit Street, 4673303  
Herzliya, Israel  
Tel: +972 9955 1344  
Fax: +972 9955 1345

**Manufacturer's Authorized  
Representative**

**Viora UK Ltd.**  
450 Bath Road, Longford  
Heathrow UB7 0EB, UK  
Tel: +44 (0) 208 757 8750  
Fax: +44 (0) 208 757 8767



This device complies with the European Medical Device Directive 93/42/EEC

V30 Platform FDA cleared for marketing  
V-FR hand piece pending FDA clearance





## Comprehensive Treatment Array

The V30 Platform houses a strong foundation of different technologies – Nd:Yag, IPL & proprietary CORE™ and SVC™ RF – that enable you to offer multiple applications to your customers. Each technology is accessible via a different handpiece, all of which can be acquired separately. You control your treatment menu.

Viora capitalizes on the advanced and proprietary technology within V30 by delivering the most effective results possible, quickly, safely, painlessly and with little to no downtime.

Below are the most popular patient needs that can be treated with the V30 Platform.



Skin Tightening

Body & Facial Contouring

Skin Rejuvenation

Skin Resurfacing

Hair Removal

Acne Clearance

Stretch Mark Reduction

Pigmented Lesions

Vascular Lesions

Cellulite Reduction

Hyperpigmentation

Reduce Fine Lines & Wrinkles



## Harnessing the Power of the V30 Platform

While utilizing single technologies with V30 results in outstanding clinically efficacious sessions, the power of the device becomes apparent when performing Viora's treatment protocols incorporating multiple technologies. Specially developed protocols founded on worldwide clinical experience allow you to deliver results at the highest standard.

## Offer More with Multiple Technologies

V30's expandable platform means that you can have the most advanced technology at your fingertips, when you want it. Interchangeable hand pieces allow you the flexibility to offer your patients a wide variety of treatments. And increasing your service offering means increasing your business potential.

### V·LASER Nd:YAG

Viora's Nd:YAG hand piece offers treatments for vascular lesions, nail fungi, skin rejuvenation and more. Four different spot sizes enable a wide variety of options, allowing you to customize treatments for each individual patient, according to his or her needs.

### V·IPL

With just one hand piece, and five interchangeable filters, you can enjoy Viora's IPL treatments, such as hair removal, skin rejuvenation, pigmentation, vascular lesions and acne clearance. Our technology ensures that every skin type and condition can be treated thanks to various pulse structures designed for safety and personalized treatments. In addition, you can further tailor your treatments with two different spot sizes for body and face.

### V·ST

The advanced ST hand piece integrates Viora's proprietary and revolutionary multi-RF technology, CORE™. The ability to do facial and body skin tightening with one hand piece is unique to Viora due to CORE's three channels of different frequency depths, and a fourth channel incorporating all three frequencies. Combined with contact cooling, Viora's ST hand piece delivers the most comfortable, effective skin tightening treatments available.

### V·FR

The second generation of fractional RF technology is solely available with the FR hand piece due to Viora's proprietary SVC™ technology. With innovative Switching, Vacuum and Cooling features, the FR hand piece makes skin resurfacing and rejuvenation a safe, effective, non-invasive procedure, with little to no downtime.



## Viora Is Invested In Your Success

Viora's commitment to you doesn't end upon your receiving your V30 Platform. We know that your success is our success and therefore we go beyond what other companies offer in clinical, marketing and technical support. Whether it's our ongoing training, advanced protocols, sales tools, marketing resources for selling to existing customers or reaching new ones, Viora puts a high priority on your needs to ensure that you receive the best and fastest return on your investment.

## Smart Business with V30

Viora's V30 is the smart acquisition for any clinic. We provide you with the most advanced technologies in a single platform so you can offer your patients the widest variety of treatments. Your possibilities with V30 are endless, as Viora is constantly developing new technologies and protocols.

# Shape the Future

